ELECTRONIC RESEARCH GROUP

P.O. Box 536 Halfway House 1685 Tel / Fax +27 11 314-1632

Email: erg@erg.co.za

Web Site: www.erg.co.za

MS2 - Mini Stage 2

OWNERS MANUAL

The **ERG Mini-Stage 2** is intended for use both in the veld and in the laboratory.

The stage surface is specially treated to enhance the heat radiation properties so that heat is transferred to the slide by both contact conduction and radiation.

This means that specimen warming is more effective when conditions are difficult such as in a cooling air flow.

In the working area the operator has the choice of using a 24mm x 6mm slot or a 15mm round hole for examination purposes.

The stage platform accommodates two slides. One on the examination area and another alongside for warming.

Temperature Control

At the end of the cast block (opposite the power cable plug) is a temperature trimming adjustment.

The stage is factory set to 37 degrees C when this control is set in the centre of the adjustment range.

The temperature can be increased or decreased by approximately 2 degrees C by turning the trimmer adjustment either clockwise or anticlockwise..

Use Of The Stage

Power Supply can be either 230v mains or 12v car battery

If the stage is to be used from a vehicle, plug the car adapter into the cigarette lighter socket and plug the stage into the adapter. A 5Amp 6x32mm fuse is inside the plug.

For use from 230 volt mains plug the mains power supply into a wall socket and plug the stage into the power supply.

The lamp on the cast block shines red when heating power is applied to the stage and green when the desired temperature has been reached.

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